

# Sequence File Syntax

- File extension - **.lseq**
- Preferred text editor – Sublime Text

## General Format

- Lines starting with **#** are considered as comments. No action will be taken on these.
- Leaving blank lines is allowed.
- Avoid any tabs or indentation in the scripts.
- All keywords must be in FULL CAPS.
- All actionable commands have the following syntax –  
`<action> <variable> <param1, param2, ...>`
- Only spaces must be used to separate words (and numbers) in actionable commands. No other punctuation should be used.

## Commands Description

- Each **action** accepts certain **variables**. *Params* list depends on the **action** and **variable** chosen.
- Each **action** is described in the following subsections.

### SET

- Drives a **variable** to a setpoint.
- Variables – **TEMP, FIELD, THETA, CURRENT, AC-SINGLE, DC-SINGLE**
- Params - `<setpoint> RATE <rate>`
- Examples –
  - **SET TEMP 300 RATE 6**
    - Set temperature to 300 **K** at a rate of 6 **K/min**
    - Do not exceed a rate of 10 K/min
  - **SET FIELD 1000 RATE 100**
    - Set field to 1000 **Oe** at a rate of 100 **Oe/s**
    - Never exceed a rate of 200 Oe/s
  - **SET THETA 90 RATE 1**
    - Set theta to 90 **degrees** at a rate of 1 **degree/s**
    - Never exceed a rate of 1 degree/s
  - **SET CURRENT 1.5e-6**
    - Set the current to  $1.5 \times 10^{-6}$  **A**
    - Note that decimal as well as exponent representation is allowed for current setpoint.
  - **SET AC-SINGLE**

- Configures instruments to perform single-shot measurements using the lock-in amplifier(s) (AC Setup settings).
- Single measurement is performed every time the **START MEASURE SINGLE** line appears in the code.
- **SET DC-SINGLE**
  - Configures instruments to perform single-shot measurements using the voltmeter(s) (DC Setup settings).
  - Single measurement is performed every time the **START MEASURE SINGLE** line appears in the code.

## WAIT

- Halts the sequence execution until setpoint of the variable is reached. Additionally delays proceeding of the sequence by specified time.
- Variables - **TEMP, FIELD, THETA, MEASURE**
- Params - *<delay time>*
- Examples –
  - **WAIT TEMP 600**
    - Waits until the temperature setpoint is reached. Then waits additional 600 **seconds** before proceeding to next lines.
    - Similar operation for **FIELD** and **THETA**
  - **WAIT MEASURE 600**
    - There is not setpoint to be reached here. This is just a provision to ensure that the sequence halts for a while until all instruments are properly initialized and ready to start meaningful measurements.
    - Halting at least 20 seconds after starting any measurement is recommended.

## SCAN – END SCAN

- Starts a loop (something like a for loop) over a variable.
- Variables - **TEMP, FIELD, THETA, CURRENT**
- Params - *<start point> <stop point> <number of points> RATE <rate>*
- Examples –
  - **SCAN FIELD 0 90000 10 RATE 100**  
**WAIT FIELD 10**  
*# Additional steps to be executed in the loop*  
**END SCAN**
    - Start a scan over field. Fields setpoints loop over 0, 10000 Oe, ..., 90000 Oe (10 points, linearly spaced, from 0 to 90000).
    - **WAIT** command right after the **SCAN** command is recommended, to allow reaching the setpoint of the current loop.
    - Do not forget the **END SCAN** command at the end of the loop.
- Note – Nested loops are NOT supported.

## START

- Start a measurement.
- Variables – `MEASURE`, `MEASURE SINGLE`, `IV`
- Params – Not applicable
- Examples –
  - `START MEASURE`
    - Starts continuous measurements based on the hardware settings (AC/DC system). Sequence execution proceeds to the next lines immediately. Measurements continue right until a `STOP MEASURE` command appears.
  - `START MEASURE SINGLE`
    - Single-shot measurement taken based on the hardware settings (AC/DC system) defined through the `SET AC-SINGLE` or `SET DC-SINGLE` commands.
  - `START IV`
    - Start an IV measurement, based on settings specified in the program. The IV plotting sub-VI pops up and closes after the IV measurement is complete. Sequence proceeds only after the IV measurement is complete.

## STOP

- Stop the measurements started with the `START MEASURE` command.
- Variables – `MEASURE`
- Params – Not applicable
- Example –
  - `STOP MEASURE`

## NEWFILE

- Create a new file for data accumulation. All data is appended in this file, over multiple `START MEASURE` – `STOP MEASURE` instances, until another file is created. File created has the chosen extension (default is `.csv`). A timestamp is also appended in the filename before the extension, in the format `_MMDD_hhmm`.
- Variables – Not applicable
- Params – *<filename without extension>*
- Examples –
  - `NEWFILE test_datafile`
    - Creates a new datafile with the name `test_datafile_1231_2359.csv` assuming `.csv` extension is chosen, and file was created on December 31 at 11:59 PM.

